

Compressive Strength Test Results

Age (Hours)	Real Compressive Strength (MPa)					
	Specimen 1	Specimen 2	Specimen 3	Specimen 4	Specimen 5	Average
1	8.86	8.01	7.42	8.07	8.14	8.10
3	14.70	15.30	15.40	14.90	14.90	15.00
7	23.80	22.50	21.90	21.70	24.10	22.80
14	29.20	29.40	28.80	27.80	29.00	28.80
28	33.20	34.10	33.60	32.60	33.60	33.40

Age (Hours)	Predicted Compressive Strength by CONCUREX					
	Specimen 6	Specimen 7	Specimen 8	Specimen 9	Specimen 10	Average
1	7.94	8.72	8.12	8.32	7.81	8.18
3	14.10	18.70	16.30	15.00	13.60	15.50
7	22.90	21.40	20.60	21.00	21.90	21.60
14	27.50	28.40	26.90	25.10	26.00	26.80
28	32.70	28.40	30.40	32.20	29.10	30.60

