

Compressive Strength Test Results

Age (Hours)	Real Compressive Strength (MPa)					
	Specimen 1	Specimen 2	Specimen 3	Specimen 4	Specimen 5	Average
0	30.20	31.30	30.10	28.00	27.90	29.50
6	66.00	70.90	66.10	71.20	69.40	68.70
12	111.00	109.00	110.00	112.00	112.00	110.80
18	138.00	139.00	141.00	142.00	139.00	139.80
24	161.00	163.00	164.00	164.00	162.00	162.80

Age (Hours)	Predicted Compressive Strength by CONCUREX (MPa)					
	Specimen 6	Specimen 7	Specimen 8	Specimen 9	Specimen 10	Average
0	27.70	28.60	27.00	27.00	33.00	28.66
6	69.80	73.40	66.90	78.90	75.00	72.80
12	108.00	109.00	110.00	101.00	98.00	105.20
18	142.00	125.00	138.00	130.00	130.00	133.00
24	157.00	148.00	165.00	151.00	166.00	157.40

